

# JOBST® ELVAREX®\*

## EFFECTIVE DAYTIME THERAPY


With a natural rubber core, Elvarex provides unsurpassed structure and control for patients of all shapes and sizes. Controlled gradient compression for the effective long-term management of lymphedema

## FLAT KNIT COMFORT

- Flat knitted construction that allows greater air permeability
- Excellent durability for extended wear



## STYLES



## COLORS

### Upper Extremity



### Lower Extremity



## COMPRESSION CLASSES

### Upper Extremity

- CCL 1 ▶ 15 - 21 mmHg
- CCL 2 ▶ 23 - 32 mmHg
- CCL 2F ▶ 23 - 32 mmHg

### Lower Extremity

- CCL 1 ▶ 18 - 21 mmHg
- CCL 2 ▶ 23 - 32 mmHg
- CCL 3 ▶ 34 - 46 mmHg
- CCL 3F ▶ 34 - 46 mmHg
- CCL 4 ▶ 49 - 70 mmHg
- CCL 4S ▶ 60 - 90 mmHg

## OPTIONS Refer to page 87

\* CAUTION: These products are made with natural rubber latex.  
 † only available in CCL1 and CCL2  
 \*\* not available in toe cap  
 Not all colors and styles may be available in all compression ranges; please refer to product charts.