

Patient Name

Measured By

Date

Clinic/DME


Purchase Order No.

Billing Address

Shipping Address

F-CWG ver.01

Comfiwave® Glove to Axilla



COMPRESSION

15-21mmHg

QUANTITY


Circumference cm

	<input checked="" type="checkbox"/> S	<input checked="" type="checkbox"/> M	<input checked="" type="checkbox"/> L
g AXILLA	25 - 30	30 - 40	40 - 50
d MID-POINT OF FOREARM	23 - 26	26 - 30	30 - 38
c WRIST	13 - 18	18 - 22	22 - 28
a PALM AT BASE FOR FINGERS	16 - 22	18 - 24	20 - 26

Length cm

	<input checked="" type="checkbox"/> SHORT	<input checked="" type="checkbox"/> REGULAR
c-g WRIST TO AXILLA	36 - 42	40 - 46

Comfiwave® Glove to Elbow



COMPRESSION

15-21mmHg

QUANTITY


Circumference cm

	<input checked="" type="checkbox"/> S	<input checked="" type="checkbox"/> M	<input checked="" type="checkbox"/> L
b BASE OF THUMB WEBBING	17 - 22	20 - 25	23 - 28
c WRIST	14 - 20	16 - 24	22 - 28
d MID - POINT OF FOREARM	23 - 28	26 - 32	30 - 38

Length cm

	<input checked="" type="checkbox"/> REGULAR
c-e WRIST TO ELBOW	20 - 24

Comfiwave® Glove



COMPRESSION

15-21mmHg

QUANTITY

Circumference cm

	<input checked="" type="checkbox"/> S	<input checked="" type="checkbox"/> M	<input checked="" type="checkbox"/> L
b BASE OF THUMB	17 - 22	20 - 25	23 - 28
c WRIST	14 - 20	16 - 24	22 - 28

