

Patient Name

Measured By

Date

Clinic/DME


Purchase Order No.

Billing Address

Shipping Address

F-EWU ver.01

### easywrap ARM



**COMPRESSION**

Light 20-30mmHg

**SIDE**

Left  Right

**COLOR**

Sand  Black

**QUANTITY**

Circumference cm

<b>g</b> AXILLA	22-31	29-39	32-45	36-50
-----------------	-------	-------	-------	-------

<b>e</b> ELBOW CREASE	20-27	25-34	30-40	32-43
-----------------------	-------	-------	-------	-------

<b>c</b> WRIST CREASE	14-18	16-21	19-25	19-25
-----------------------	-------	-------	-------	-------


**S**
 **M**
 **L**
 **XL**

Length cm

<b>c-g</b> WRIST TO AXILLA	40-44	44-48	48-52
----------------------------	-------	-------	-------

**SHORT**
 **REGULAR**
 **LONG**

### easywrap HAND



**COMPRESSION**

Light 20-30mmHg

**SIDE**

Left  Right

**COLOR**

Sand  Black

**QUANTITY**

Circumference cm

<b>b</b> BASE OF THUMB	<19	19-21	21-23	>24
------------------------	-----	-------	-------	-----

<b>a</b> PALM AT BASE OF FINGERS	<18	18-20	20-22	>22
----------------------------------	-----	-------	-------	-----

**S**
 **M**
 **L**
 **XL**



## Left

### easywrap<sup>®</sup> ARM

COLOR	LENGTH	S	M	L	XL
SAND	SHORT	EL-AS2L	EL-AS3L	EL-AS4L	EL-AS5L
	REGULAR	EL-AR2L	EL-AR3L	EL-AR4L	EL-AR5L
	LONG	EL-AL2L	EL-AL3L	EL-AL4L	EL-AL5L
BLACK	SHORT	ELB-AS2L	ELB-AS3L	ELB-AS4L	ELB-AS5L
	REGULAR	ELB-AR2L	ELB-AR3L	ELB-AR4L	ELB-AR5L
	LONG	ELB-AL2L	ELB-AL3L	ELB-AL4L	ELB-AL5L

### easywrap<sup>®</sup> HAND

COLOR	XS	S	M	L
SAND	EL-H1L	EL-H2L	EL-H3L	EL-H4L
BLACK	ELB-H1L	ELB-H2L	ELB-H3L	ELB-H4L

## Right

S	M	L	XL
EL-AS2R	EL-AS3R	EL-AS4R	EL-AS5R
EL-AR2R	EL-AR3R	EL-AR4R	EL-AR5R
EL-AL2R	EL-AL3R	EL-AL4R	EL-AL5R
ELB-AS2R	ELB-AS3R	ELB-AS4R	ELB-AS5R
ELB-AR2R	ELB-AR3R	ELB-AR4R	ELB-AR5R
ELB-AL2R	ELB-AL3R	ELB-AL4R	ELB-AL5R

XS	S	M	L
EL-H1R	EL-H2R	EL-H3R	EL-H4R
ELB-H1R	ELB-H2R	ELB-H3R	ELB-H4R

