

JOBST® CONFIDENCE

NEXT GENERATION OF MADE-TO-MEASURE FLAT-KNIT COMPRESSION GARMENTS

- 3D innovative knitting method allows to conform to individual body shapes
- Conforming fabric that also offers firm support
- Moisture wicking with double-layer construction for optimal comfort
- Flexible fabric to make putting on and taking off easier
- Not made with natural rubber latex



STYLES



COLORS



COMPRESSION CLASSES

Upper Extremity

- CCL 1 ▶ 15 - 21 mmHg
- CCL 2 ▶ 23 - 32 mmHg

Lower Extremity

- CCL 1 ▶ 18 - 21 mmHg
- CCL 2 ▶ 23 - 32 mmHg
- CCL 3 ▶ 34 - 46 mmHg

JOBST® Confidence Design Options on next page